



## **IA2 ELISA**

For the quantitative and qualitative determination of antibodies  
against human islet antigen 2 (IA2) in human serum

For Research Use Only. Not For Use In Diagnostic Procedures.

**Catalog Number:** 35-IA2HU-E01

**Size:** 96 wells

**Version:** RUO: 2023-03-02 (V004) – ALPCO 1.0

## 1. Intended Use

The IA2 ELISA is a solid phase enzyme immunoassay employing recombinant human tyrosine phosphatase-related protein islet antigen 2 (IA2) for the quantitative and qualitative detection of antibodies against human IA2 in human serum. For Research Use Only. Not for use in Diagnostic Procedures.

## 2. Principle of the Assay

Serum samples are diluted 1:4 and incubated in the microtiter plates coated with the specific antigen. Antibodies, if present in the sample, bind to the antigen. The unbound fraction is washed off in the following step. Afterwards biotin-labeled IA2 is incubated. During this step IA2 autoantibodies form a bridge between the IA-2 immobilized on the plate and the IA2-Biotin in the liquid phase. Unbound IA2 is washed off in the following step. The amount of bound IA2-Biotin is then determined by the addition of streptavidin-peroxidase (conjugate) that binds to biotin. Unbound streptavidin-peroxidase is washed off in the following step. Addition of TMB-substrate generates an enzymatic colorimetric (blue) reaction, which is stopped by diluted acid (color changes to yellow). The rate of color formation from the chromogen is proportional to the initial concentration of the respective antibodies in the sample.

## 3. Kit Contents

<b>To Be Reconstituted</b>				
<b>Item</b>	<b>Quantity</b>	<b>Cap Color</b>	<b>Solution color</b>	<b>Description/Contents</b>
Sample Buffer (5x)	1 x 20ml	White	Yellow	5x concentrated Tris, sodium chloride (NaCl), bovine serum albumin (BSA), sodium azide < 0.1% (preservative)
Wash Buffer (50x)	1 x 20ml	White	Green	50x concentrated Tris, NaCl, Tween 20, sodium azide < 0.1% (preservative)
IA2 Biotin	3 vials	White	-	Biotin conjugated to recombinant human IA2; bovine serum albumin (BSA), lyophilized
<b>Ready to Use</b>				
<b>Item</b>	<b>Quantity</b>	<b>Cap color</b>	<b>Solution color</b>	<b>Description/Contents</b>
Reconstitution Buffer	1 x 20ml	Red	Red	PBS, bovine serum albumin (BSA)
Negative Control	1 x 1ml	Green	Colorless	Control material (diluted), bovine serum albumin (BSA), sodium azide < 0.1% (preservative)
Positive Control	1 x 1ml	Red	Yellow	Control material (diluted), bovine serum albumin (BSA), sodium azide < 0.1% (preservative)
Cut-off Calibrator	1 x 1ml	Blue	Yellow	Calibrator material (diluted), bovine serum albumin (BSA), sodium azide < 0.1% (preservative)

Calibrators	6 x 1ml	White	Yellow *	Concentration of each calibrator: 0, 25, 75, 125, 250, 500 IU/ml. Calibrator material (diluted), bovine serum albumin (BSA), sodium azide < 0.1% (preservative)
Conjugate	1 x 15ml	Blue	Blue	Streptavidin coupled to horseradish peroxidase; bovine serum albumin (BSA)
TMB Substrate	1 x 15ml	Black	Colorless	Stabilized tetramethylbenzidine and hydrogen peroxide (TMB/H <sub>2</sub> O <sub>2</sub> )
Stop Solution	1 x 15ml	White	Colorless	1M Hydrochloric Acid
Microtiter plate	12 x 8 well strips	N/A	N/A	With breakaway microwells. Refer to paragraph 1 for coating.
• Color increasing with concentration				
<b>Materials Required, but not provided</b>				
Microtiter plate reader 450 nm reading filter and recommended 620 nm reference filter (600-690 nm). Glassware (cylinder 100-1000ml), test tubes for dilutions. Vortex mixer, precision pipettes (10, 100, 200, 500, 1000 µl) or adjustable multipipette (100-1000µl). Microplate washing device (300 µl repeating or multichannel pipette or automated system), adsorbent paper. Our tests are designed to be used with purified water according to the definition of the United States Pharmacopeia (USP 26 - NF 21) and the European Pharmacopeia (Eur.Ph. 4th ed.).				

#### 4. Storage and Shelf Life

Store all reagents and the microplate at 2-8°C/35-46°F, in their original containers. Please note: Reconstituted IA2 Biotin must be used on day of reconstitution. Once prepared, other reconstituted solutions are stable at 2-8°C/35-46°F for 1 month. Reagents and the microplate must be used within the expiration date indicated on each component. Avoid exposing TMB solution to light. Store microplates in designated foil, with desiccant, and seal tightly.

#### 5. Precautions

##### 5.1 Health hazard data

This product is for research use only. Not for use in diagnostic procedures. Only staff trained and specially advised in the methods of ELISA techniques may perform the kit. Although this product is not considered particularly toxic or dangerous in conditions of normal use, refer to the following for maximum safety:

##### **Recommendations and precautions**

This kit contains potentially hazardous components. Though kit reagents are not classified as being irritants to eyes and skin it is recommended to avoid contact with eyes and skin and to wear disposable gloves.

**WARNING!** Calibrators, controls, and buffers contain sodium azide (NaN<sub>3</sub>) as a preservative. NaN<sub>3</sub> may be toxic if ingested or absorbed by skin or eyes. NaN<sub>3</sub> may react with lead and copper plumbing to form highly explosive metal azides. Upon disposal, flush with plenty of water to prevent azide build-up. Please refer to decontamination procedures as outlined by the CDC or other local/national guidelines.

**Do not smoke, eat, or drink when manipulating the kit. Do not pipette by mouth.**

All biological source material used in kit reagents has been tested by approved methods and found negative for HBsAg, Hepatitis C, and HIV-1. However, no test can guarantee the absence of viral agents in such material completely. Handle kit controls, standards, and samples as if capable of transmitting infectious diseases and according to local and national requirements.

This kit contains material of animal origin as stated in the materials included, handle according to local and national requirements.

## 5.2 General directions for use

- In case that the product information, including the labeling, is defective or incorrect please contact ALPCO.
- Do not mix or substitute reagents or microplates from different lot numbers. This may lead to variations in the results.
- Allow all components to reach room temperature (20-32°C/68-89.6°F) before use, mix well, and follow the recommended incubation scheme for optimum performance of the assay.
- **Incubation: It is recommended to test performance at 30°C/86°F for automated systems.**
- Never expose components to higher temperatures than 37°C/98.6 °F.
- Always pipette substrate solution with brand new tips only. Protect this reagent from light. Never pipette conjugate with tips used previously with other reagents.

## 6. Sample Collection, Handling, and Storage

It is preferable to use freshly collected serum samples. Blood draw must follow local and national requirements. Do not use icteric, lipemic, hemolyzed, or bacterially contaminated samples. Sera with particles should be cleared by low-speed centrifugation (<1000 x g). Blood samples should be collected in clean, dry, and empty tubes.

After separation, the serum samples should be used during the first 8 hours, respectively stored tightly closed at 2-8°C/35-46°F for up to two days (48 hours), or frozen at -20°C/-4°F for longer periods. (Thomas: Labor und Diagnose; CLSI Guideline GP44-A4).

## 7. Assay Procedure

### 7.1 Preparations prior to starting

#### Dilute concentrated reagents:

- Dilute the concentrated sample buffer 1:5 with distilled water (e.g., 20 mL plus 80 mL).
- Dilute the concentrated wash buffer 1:50 with distilled water (e.g., 20 mL plus 980 mL).
- To avoid mistakes, it is suggested to mark the cap of the different calibrators.

#### Samples:

Dilute serum samples 1:4 with 1x working sample buffer e.g. 225 µl sample buffer (1x) + 75 µl serum. Mix well!

#### Reconstitute and dilute IA2 Biotin:

Reconstitute lyophilized IA2 Biotin in two steps: First, add 1 ml reconstitution buffer to one vial, let it stand for 5 minutes and mix it until IA2 Biotin is completely dissolved. Make sure that no remaining lyophilized or reconstituted IA2 Biotin is in the cap. In a second step, transfer 1 ml reconstituted IA2 Biotin to 5 ml reconstitution buffer and mix well. One vial of reconstituted IA2 Biotin is sufficient for six strips. If more strips will be used in one test, the two vials of lyophilized IA2 Biotin must be reconstituted and diluted as described above and then combined. Mix well!

**Washing:**

Prepare 20 mL of 1x working wash buffer per 8 wells or 200 mL for 96 wells

- e.g., 4 mL concentrate plus 196 mL distilled water.

**Automated washing:**

Consider excess volumes required for setting up the instrument and the dead volume of robotic pipette.

**Manual washing:**

Discard liquid from wells by inverting the plate. Knock the microwell frame with wells downside vigorously on clean absorbent paper. Pipette 300 µL of 1x working wash buffer into each well, wait for 20 seconds. Repeat the whole procedure two more times.

**Microplates:**

Calculate the number of wells required for the test. Remove unused wells from the frame. Replace and store in the provided plastic bag, together with desiccant, seal tightly (2-8°C/35-46°F).

**7.2 Suggested Layout**

The following layout for pipetting the calibrators, controls, and samples is suggested:

	For quantitative interpretation				For qualitative interpretation			
	1	2	3	4...	1	2	3	4...
<b>A</b>	Cal A	Cal E	S1		NC	S2		
<b>B</b>	Cal A	Cal E	S1		NC	S2		
<b>C</b>	Cal B	Cal F	S2		CC	S3		
<b>D</b>	Cal B	Cal F	S2		CC	S3		
<b>E</b>	Cal C	PC	S3		PC	...		
<b>F</b>	Cal C	PC	S3		PC			
<b>G</b>	Cal D	NC	...		S1			
<b>H</b>	Cal D	NC	...		S1			

Cal A: calibrator A      Cal D: calibrator D      PC: Positive Control      S1: Sample 1  
 Cal B: calibrator B      Cal E: calibrator E      NC: Negative Control      S2: Sample 2  
 Cal C: calibrator C      Cal F: calibrator F      CC: Cut-Off Calibrator      S3: Sample 2

**7.3 Assay Steps**

1. Ensure preparations from step 7.1 have been completed prior to pipetting.
2. Use the following steps in accordance with the quantitative/qualitative results desired:
3. Pipette, into the designated microwells as described in section 7.2, 100 µL of either:
  - a. Calibrators (Cal A to Cal F) for Quantitative or
  - b. Cut-off Calibrator (CC) for Qualitative

And 100 µL of each of the following:

- c. Negative control (NC) and Positive Control (PC) and
- d. Diluted serum samples (S1, S2...)
4. Incubate for 60 minutes at 20-32°C/68-89.6°F.
5. Wash 3 times with 300 µL of 1x working wash buffer.
6. Pipette 100 µL of IA2 Biotin solution (as described in chapter 7.1 above) into each well.
7. Incubate for 30 minutes at 20-32°C/68-89.6°F.
8. Wash 3 times with 300 µL of 1x working wash buffer.
9. Pipette 100 µL of conjugate into each well.
10. Incubate for 30 minutes at 20-32°C/68-89.6°F.
11. Wash 3 times with 300 µL of 1x working wash buffer.
12. Pipette 100 µl of TMB Sustrate into each well.
13. Incubate for 30 minutes at 20-32°C/68-89.6°F, protected from light.
14. Pipette 100 µl of stop solution into each well, using the same order as when the substrate was pipetted.
15. Incubate plate for a minimum of 5 minutes.
16. Agitate plate carefully for 5 seconds by tapping the plate frame.
17. Read absorbance at 450 nm (450/620 nm if available) within 30 minutes.

## 8. Calculation of Quantitative and Qualitative Results

For **quantitative interpretation** establish the standard curve by plotting the **optical density (OD) of each standard (y-axis)** with respect to the corresponding concentration values in **IU/mL (x-axis)**. Log/lin coordinates and 4-Parameter Fit are recommended for best results. From the OD of each sample, read the corresponding antibody concentrations expressed in **IU/mL**.

Normal Range	Cut-off Range	Positive Range
< 20 IU/mL	20 - 30 IU/mL	>30 IU/mL

### **Example of a standard curve**

Do not use this data for calculating assay results.

Calibrators IgG	OD 450/620 nm	CV % (Variation)
0 IU/ml	0.072	1.9
25 IU/ml	0.241	1.8
75 IU/ml	0.599	2.7
125 IU/ml	1.017	3.6
250 IU/ml	1.382	5.0
500 IU/ml	2.145	1.2

### Example calculation

Sample	Replicate (OD)	Mean (OD)	Result (IU/mL)
S 1	0.808/0.831	0.820	103.5
S 2	1.081/1.071	1.076	148.1

- Samples above the highest calibrator range should be reported as > Max and diluted as appropriate and re-assayed. Samples below the calibrator range should be reported as < Min.
- For lot-specific data, see enclosed quality control leaflet. Laboratories might perform an in-house quality control by using their own controls and/or internal pooled sera, as foreseen by national regulations.
- Each laboratory should establish its own normal range according to its own established procedures.
- If the values of the controls do not meet the criteria, the test is invalid and must be repeated.
- The following technical issues should be verified: Expiration dates of (prepared) reagents, storage conditions, pipettes, devices, photometer, incubation conditions, and washing methods.
- If the items tested show aberrant values or any kind of deviation or that the validation criteria are not met without explicable cause, please contact ALPCO.
- **This kit is for research use only. It is not for use in diagnostic procedures.**

For qualitative interpretation read the optical density of the cut-off control and the samples. Compare the sample's OD with the OD of the cut-off calibrator. For qualitative interpretation, it is recommended to consider sera within range of 20% around the cut-off value as equivocal. All samples with higher ODs are considered positive, samples with lower ODs are considered negative.

Negative:	$OD_{\text{sample}} < 0.8 \times OD_{\text{cut-off}}$	
Equivocal:	$0.8 \times OD_{\text{cut-off}} \leq OD_{\text{sample}} \leq 1.2 \times OD_{\text{cut-off}}$	
Positive:	$OD_{\text{sample}} > 1.2 \times OD_{\text{cut-off}}$	

## 9. Assay Characteristics

### 9.1 Summary

Sample material:	serum
Sample volume:	75 µl of sample diluted 1:4 with 1x sample buffer
Total incubation time:	150 minutes at 20-32°C/68-89.6°F
Calibration range:	0 - 500 IU/ml
Analytical sensitivity:	3.98 IU/ml
Storage:	at 2-8°C/35-46°F in original vials only
Number of determinations:	96 tests

## 10. Literature

1. **Chen S, Willis J, Maclean C, Ananieva-Jordanova R, Amoroso MA, Brooking H, Powell M, Collins A, Bennett S, Mitchell S, Burne P, Furmaniak J, Smith BR (2005).** Sensitive non-isotopic assays for autoantibodies to IA-2 and to a combination of both IA-2 and GAD65. *Clin Chim Acta.*357 (1): 74-83.
2. **Törn C, Mueller PW, Schlosser M, Bonifacio E, Bingley PJ; Participating Laboratories (2008).** Diabetes Antibody Standardization Program: evaluation of assays for autoantibodies to glutamic acid decarboxylase and islet antigen-2. *Diabetologia* 51 (5): 846-852.
3. **Lan MS, Wasserfall C, Maclaren NK, Notkins AL (1996).** IA-2, a transmembrane protein of the protein tyrosine phosphatase family, is a major autoantigen in insulin-dependent diabetes mellitus. *Proc Natl Acad Sci U S A.*; 93 (13): 6367-6370.
4. **Van den Driessche A, Eenkhoorn V, Van Gaal L, De Block C (2009).** Type 1 diabetes and autoimmune polyglandular syndrome: a clinical review. *Neth J Med.* 67(11):376-387.
5. **Winter WE, Harris N, Schatz D (2002).** Type 1 diabetes islet autoantibody markers. *Diabetes Technol Ther* 4(6):817-839.
6. **Winter WE (2015).** IASP 2015 - The current workshop and planning of future workshops. Pre-sented at the 14th International Congress of the Immunology of Diabetes Society (12-16 April 2015, Munich, Germany).
7. **Lampasona V, Williams AJK, Schlosser M (2015).** Results of the IASP 2015 Workshop. Presented at the 14th International Congress of the Immunology of Diabetes Society (12-16 April 2015, Munich, Germany).
8. **Lothar Thomas:** Labor und Diagnose. Indikation und Bewertung von Laborbefunden für die medizinische Diagnostik., 8. Auflage, TH Books
9. **CLSI Guideline GP44-A4:** Procedures for the Handling and Processing of Blood Specimens for Common Laboratory Tests